

AutoTrace 280: Solid-Phase Extraction Instrument



The AutoTrace® 280 instrument handles large volume Solid-Phase Extraction (SPE) of liquid samples from 20 mL to 4 liters. The instrument can be configured with 1, 3, or 6 mL SPE cartridges or 47 mm solid-phase disks. The instrument automatically loads liquid matrices on the SPE material then automatically elutes the analytes of interests using up to five different solvents.

System Features

- Automatically conditions, rinses, and elutes SPE cartridges with a choice of five solvents or reagents
- Simultaneously loads up to six samples
- Loads sample volumes from 20 mL to 4000 mL
- Offers choice of four different collection vial racks
- Segregates aqueous and solvent waste using separate waste lines
- Processes ground, surface, and wastewater samples – even those containing suspended solids
- Easy-to-use SPE specific software
- Supports SPE cartridges or discs
- Documented EPA methods for water testing
- Integrated venting of solvent vapors

Applications

The AutoTrace instrument is perfect for laboratories looking to automate large volume water extractions. AutoTrace applications include many of the U.S. EPA methods and a variety of other SPE sample preparation cleanup methods.

Methods

- U.S. EPA Method 525
- U.S. EPA Method 549.2
- U.S. EPA Method 508 & 608
- U.S. EPA Method 625
- U.S. EPA Method 1664
- Extraction of total petroleum hydrocarbons
- Extraction of acid herbicides
- Extraction of triclosan
- Extraction of nitrosamines (EPA Method 521)
- Extraction of steroids and phytoestrogens
- Extraction of endocrine disruptors (EPA Method 1694)
- Extraction of personal care products



Passion. Power. Productivity.

TABLE 1. PERFORMANCE				
Pesticide Recovery Study N = 6	AutoTrace SPE		Vacuum Manifold SPE	
	Recovery%	RSD	Recovery%	RSD
Atrazine	88	1.8	54	12.2
Propazine	91	1.5	80	7.3
Alachor	99	3.4	96	4.1
Metachor	99	4.3	96	2.9

Table 1 shows a pesticide recovery study comparing the AutoTrace SPE technique and a vacuum manifold technique. The improvements in recovery and reproducibility are attributed to the microprocessor control of all the liquid flow rates – both sample and SPE reagents.

KEY SPECIFICATIONS	
<p>Gas Regulator and Gas Gauge Range: Output: 0–30 psi (0–1.4 bars), Input: 100 psi (6.9 bars) maximum</p> <p>Net Weight: 150 lbs. (68.1 kg)</p> <p>Dimensions (h × w × d): 23 × 25 × 27 in (57 × 63.5 × 69 cm) Sample Rack: 8 × 16.5 × 13 in)</p> <p>500 & 600 Series EPA Methods: 24 methods stored in the AutoTrace instrument, unlimited number of methods stored in PC</p> <p>Electrical: Voltage: 100, 120, 220, or 240 V ± 10% Frequency: 47–63 Hz Power: 150 Volt AC</p> <p>Liquid Management: Air Syringe: One 10 mL air syringe LH Syringe: One 10 mL liquid handling syringe 12-port Valve: Rotary, sliding Rulon® seal Valves: 3-way, Teflon® Nozzles: Stainless steel Tubing: TFE, 1/16" ID, (6 total supplied)</p>	<p>Sample Pumps: Displacement: Positive Accuracy: ± 2.5% Tube Fitting: Kynar® Piston and Liner: Ceramic Non Use: Not for use with acetic acid and acetone</p> <p>Solid Phase Extraction Configurations: 1 mL Syringe: Compatible cartridge 3 mL Syringe: Compatible cartridge 6 mL Syringe: Compatible cartridge 47 mm SPE disk</p>

ORDERING INFORMATION	
To order, use the following part numbers and contact your local Dionex office or distributor nearest you. In the U.S., call (800) 346-6390. In other regions, refer to the phone numbers below.	
Description	
AutoTrace 280 Automated Large Volume SPE for Disks	071386
AutoTrace 280 Automated Large Volume SPE for 6 mL Cartridges	071385
AutoTrace 280 Automated Large Volume SPE for 3 mL Cartridges	072605
AutoTrace 280 Automated Large Volume SPE for 1 mL Cartridges	072604
AutoTrace 280 Automated Large Volume SPE for 6 mL Glass Cartridge.....	072606

AutoTrace is a registered trademark of Dionex Corporation.
Rulon is a registered trademark of Saint Gobain. Teflon is a registered trademark of E.I. du Pont de Nemours and Company.
Kynar is a registered trademark of Arkema Inc.



Dionex Corporation
1228 Titan Way
P.O. Box 3603
Sunnyvale, CA
94088-3603
(408) 737-0700

North America
U.S./Canada (847) 295-7500

South America
Brazil (55) 11 3731 5140

Europe
Austria (43) 1 616 51 25 Benelux (31) 20 683 9768; (32) 3 353 4294
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

Asia Pacific
Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190
Taiwan (886) 2 8751 6655

www.dionex.com



LPN 2160-01 5M 08/09
©2009 Dionex Corporation