

PSX 40

Features

- excellent price-performance ratio
- use of two well plates (up to 384 wells) or plastic sample racks (variable, up to 50 ml)
- 15 different rack versions are applicable (codable)
- typically sample vials 2 ml, 4 ml, 8 ml
- integrated wash position (different solutions)
- integrated peristaltic pump
- 3 external positions
- diluter can be connected by RS232
- display for status information (i.e. current sample position)
- control via PC or via your analytical system
- individual system configurations are possible
- also available as OEM-modul

Technical data:

- USB-interface
- 24 V DC (external power supply 110-230 V AC)
- w x h x d: 310 x 215 x 200 mm
- weight: 4 kg
- port for diluter
- port for 8-Port-Distribution valve or injection valve



MLE

flexible XY-Autosampler

small and easy in handling

applicable in biotechnological and chemical laboratories

- automated sample handling
- use as an autosampler or sample collector
- supports different measurement routines
- easy adaption to your analytical system
- dilution function
- 384 well plates usable
- minimal space requirements



Measurement

Laboratory

Engineering