

# SofTop™ Quik Dispenser

## The Ergonomic Solution in Reagent Dispensing

- Safe, simple operation with ergonomic sure-grip design
- Maximum chemical and thermal resistance with innovative ceramic plunger
- Fast and accurate micrometer volume adjustment

The SofTop Quik Dispenser offers an exceptional “soft feel” for precision liquid handling from the creator of the world’s first ergonomic pipette. The graduated gauge offers user independent reproducible volumetric delivery. Calibrated volumetric steps assure identical volume delivery every time.

The ceramic plunger is manufactured from aluminum oxide ceramic for dimensional stability and functional safety. With a ceramic plunger, the chance of abrasions, scratching, swelling or flaking of a coating layer is completely eliminated.

The SofTop Quik Dispenser has a screw-fitted suction hose, eliminating the possibility of air bubbles during aspiration. The fully autoclavable SofTop Quik Dispenser is available in six different volumes from 0.2 mL to 60 mL. Each dispenser is individually tested and is sold with a Certificate of Quality.



## Technical Data

Volumes, mL	Increments, mL	Accuracy, %	Reproducibility, %
0.2 – 1.0	0.05	0.6	0.1
0.4 – 2.0	0.05	0.5	0.1
1.0 – 5.0	0.10	0.5	0.1
2.0 – 10.0	0.25	0.5	0.1
5.0 – 30.0	0.50	0.5	0.1
10.0 – 60.0	1.00	0.5	0.1

Typical dispensing temperatures: 4°C to 50°C.

## Ordering Information

Description	Part #
SofTop Quik Dispenser 0.2 – 1.0 mL	51249-01
SofTop Quik Dispenser 0.4 – 2.0 mL	51249-20
SofTop Quik Dispenser 1.0 – 5.0 mL	51249-50
SofTop Quik Dispenser 2.0 – 10 mL	51249-10
SofTop Quik Dispenser 5.0 – 30 mL	51249-30
SofTop Quik Dispenser 10 – 60 mL	51249-60

## Replacement Parts for SofTop™ Quik Dispensers

Description	Part #
Dispensing Nozzle for SofTop Quik Dispensers up to 10 mL	51249-70
Dispensing Nozzle for SofTop Quik Dispensers 30 and 60 mL	51249-72
12" Suction Tube for SofTop Quik Dispensers up to 10 mL	51249-74
12" Suction Tube for SofTop Quik Dispensers 30 and 60 mL	51249-76

Call us today to learn more about the SofTop™ Quik Dispenser, the newest member of Hamilton's family of ergonomic products that are easy-to-use, reliable, and maximize your laboratory's efficiency.



**HAMILTON**  
THE MEASURE OF EXCELLENCE™

Hamilton Company  
4970 Energy Way  
Reno, Nevada 89502 USA  
Toll-Free: 800-648-5950  
Telephone: +1-775-858-3000  
Fax: +1-775-856-7259  
e-mail: sales@hamiltoncompany.com

<http://www.hamiltoncompany.com>

Lit. No. L70024 ©Hamilton Company 1/05 5K Printed in U.S.A.