


| | E-1000V | E-1000UV |
|---|--|--|
| Model |  | |
| Display | 70*40mm blue backlit LCD | |
| Wavelength Range | 320 - 1020nm | 190 - 1020nm |
| Slit Width | 4 nm | |
| Wavelength Accuracy | ±2nm | |
| Wavelength Repeatability | ≤1nm | |
| Photometric Accuracy | 0.5%T | |
| Photometric Repeatability | 0.2%T | |
| Stray light | ≤0.15%T@ 360nm | |
| Stability | 0.002A@ 500nm | |
| Output Port | RS232 | |
| Light Source | Tungsten Halogen Lamp | Tungsten Halogen/Deuterium Lamp |
| Power Requirements | 110-220V/50-60HZ | |
| Shipping Dimensions & Weight | 520*450*300mm ,9KG | |