

GC Column Cross-Reference: Columns by Phase

SCION	Phase Description	Phenomenex	Agilent	Alltech	Supelco	SGE	Restek
SCION-1	dimethyl polysiloxane	ZB-1	HP-1, DB-1, CP-Sil 5 CB	AT-1, EC-1	SPB-1	BP1	Rtx-1
SCION-1MS	dimethyl polysiloxane (low bleed)	ZB-1, ZB-1ms	HP-1ms, HP-1ms UI, DB-1ms, DB-1ms UI, VF-1ms, Ultra-1	AT-1ms	SPB-1, Equity-1	BP1	Rxi-1ms
SCION-5	diphenyl dimethyl polysiloxane	ZB-5	HP-5, DB-5, CP-Sil 8 CB	EC-5, AT-5	SPB-5	BP5	Rtx-5 Rtx-5MS
SCION-5MS	diphenyl dimethyl polysiloxane (low bleed)	ZB-5, ZB-5msi	HP-5ms, HP-5ms UI, DB-5, Ultra-2, CP-Sil 8 CB	AT-5ms	SPB-5, Equity-5	BP5ms	Rxi-5ms
SCION-35	diphenyl dimethyl polysiloxane	ZB-35	HP-35, DB-35	AT-35, AT-35-ms	SPB-35, SPB-608		Rtx-35
SCION-35MS	diphenyl dimethyl polysiloxane	MR2	DB-35ms, DB35msUI, VF-35ms			BPX35, BPX608	Rxi-35Sil MS
SCION-624MS	cyanopropylmethyl phenylmethyl polysiloxane	ZB-624	DB-1301, DB-624, DB-624UI, VF-1301ms, VF-624ms, CP-1301	AT-624, AT-1301	SPB-624	BP624	Rtx-1301 Rtx-624
SCION-1701	cyanopropylmethyl phenylmethyl polysiloxane	ZB-1701, ZB-1701P	DB-1701P, DB-1701, CP-Sil 19 CB, VF-1701ms, VF-1701 Pesticides	AT-1701	Equity-1701	BP10	Rtx-1701
SCION-WAX	polyethylene glycol	ZB-Wax	DB-Wax	AT-WAX, EC-WAX		BP20	Rtx-WAX Stabilwax
SCION-WAXMS	polyethylene glycol	ZB-WAXplus	HP-INNOWax, CP-Wax 52 CB, VF-WAX MS	AT-WAXms, EC-Wax	Supelcowax- 10		Stabilwax-MS

SCION Instruments

1 Bain Square
Kirkton Campus
Livingston, Scotland
EH54 7DQ
United Kingdom
+44 (0) 1506 300 200
sales-eu@scioninstruments.com

SCION Instruments NL BV
Stanleyweg 4
Goes
4462GN
Netherlands
+ 31 (0) 113 348926
sales-eu@scioninstruments.com